

## KELLY RETRACTOR 23.5 CM - 40X50 MM

Réf. IE02701



## KELLY RETRACTOR 23.5 CM - 40X50 MM



Réf. IE02701

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040084619
Guarantee (month):	24	Customs Code	9018908500
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	3,2	Width packaging (cm):	
Length (cm) :	25,5	Length packaging (cm) :	
Height (cm) :	5,8	Height packaging (cm) :	
Weight (kg) :	0,108	Weight packaging (kg) :	0,108
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



19/05/2022 21:10 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IE02701

All characteristics are given for information only and may be changed