

HARTMANN RETRACTOR STAINLESS STEEL WIRE 16 CM (PAIR)

Réf. IE03616



HARTMANN RETRACTOR STAINLESS STEEL WIRE 16 CM (PAIR)



Réf. **IE03616**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040016375
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	5,5	Width packaging (cm):	
Length (cm) :	16,0	Length packaging (cm) :	
Height (cm) :	4,0	Height packaging (cm) :	
Weight (kg) :	0,085	Weight packaging (kg) :	0,100
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:30 - This article is available on our website : www.gsh-med.fr with the following reference IE03616

All characteristics are given for information only and may be changed