

## CUSHING RETRACTOR 20 CM 15 MM

Réf. IE08820



## CUSHING RETRACTOR 20 CM 15 MM



Réf. IE08820

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040016429
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,4	Width packaging (cm):	
Length (cm) :	20,0	Length packaging (cm) :	
Height (cm) :	1,0	Height packaging (cm) :	
Weight (kg) :	0,040	Weight packaging (kg) :	0,040
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:31 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IE08820

All characteristics are given for information only and may be changed