

DOYEN VAGINAL SPECULUM 90 MM

Réf. IV02690



DOYEN VAGINAL SPECULUM 90 MM



Réf. IV02690

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040032306
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	14,5	Width packaging (cm):	
Length (cm) :	25,2	Length packaging (cm) :	
Height (cm) :	4,6	Height packaging (cm) :	
Weight (kg) :	0,108	Weight packaging (kg) :	0,109
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:50 - This article is available on our website : www.gsh-med.fr with the following reference IV02690

All characteristics are given for information only and may be changed