

DOYEN VAGINAL SPECULUM 120 MM

Réf. IV02720



DOYEN VAGINAL SPECULUM 120 MM



Réf. IV02720

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040032313
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	12,3	Width packaging (cm):	
Length (cm) :	25,0	Length packaging (cm) :	
Height (cm) :	4,8	Height packaging (cm) :	
Weight (kg) :	0,216	Weight packaging (kg) :	0,216
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:50 - This article is available on our website : www.gsh-med.fr with the following reference IV02720

All characteristics are given for information only and may be changed