

KIRSCHNER VALVE 90° 16 CM - 105 MM X 38 MM

Réf. IV06686



KIRSCHNER VALVE 90° 16 CM - 105 MM X 38 MM



Réf. IV06686

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040134215
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,151	Weight packaging (kg) :	0,151
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:27 - This article is available on our website : www.gsh-med.fr with the following reference IV06686