

ELITE BAGS ISOTHERM BAG/CASE - COOL

ELITE BAGS

Réf. MECOOL03_B



GARANTIE
2
ANS
YEARS
GARANTEE



ELITE BAGS ISOTHERM BAG/CASE - COOL



Réf. MECOOL03_B

FEATURES			
Medical device	No		
Marquage	Barcode	8436564540355	
Guarantee (month):	24	Customs Code	4202.12.99.90
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	29,0	Width packaging (cm):	30,0
Length (cm) :	44,0	Length packaging (cm) :	44,5
Height (cm) :	39,0	Height packaging (cm) :	38,0
Weight (kg) :	5,300	Weight packaging (kg) :	5,100
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:18 - This article is available on our website : www.gsh-med.fr with the following reference MECOOL03_B