

## ELITE BAGS-EXTREME-POLYESTER/WATERPROOF- RED


Réf. MEEXTR02\_R

 ELITE BAGS

GARANTIE  
**2**  
ANS  
YEARS  
GUARANTEE



CE



## ELITE BAGS-EXTREME-POLYESTER/WATERPROOF- RED



Réf. MEEEXTR02\_R

FEATURES			
Medical device	No		
Marquage	CE	Barcode	8436564540300
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	22,0	Width packaging (cm):	23,5
Length (cm) :	47,0	Length packaging (cm) :	43,0
Height (cm) :	28,0	Height packaging (cm) :	39,0
Weight (kg) :	2,860	Weight packaging (kg) :	3,800
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



12/04/2022 21:08 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference MEEEXTR02\_R

2/2