


## GALEZOWSKI LACHRYMAL PROBE BULBOUS TIP 10 CM N°6 DIA 1.5 MM

Réf. IS05106



# GALEZOWSKI LACHRYMAL PROBE BULBOUS TIP 10 CM N°6 DIA 1.5 MM



Réf. IS05106

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040030388
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,2	Width packaging (cm):	
Length (cm) :	11,2	Length packaging (cm) :	
Height (cm) :	0,1	Height packaging (cm) :	
Weight (kg) :	0,004	Weight packaging (kg) :	0,004
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:06 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IS05106

All characteristics are given for information only and may be changed