

FLIP COVER DORO 8035, BLACK



Réf. HDACCE128N



FLIP COVER DORO 8035, BLACK



Réf. HDACCE128N

FEATURES			
Medical device	No		
Marquage		Barcode	7322460074586
Guarantee (month):	0	Customs Code	4202321000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	11,5
Length (cm) :		Length packaging (cm) :	19,0
Height (cm) :		Height packaging (cm) :	9,3
Weight (kg) :		Weight packaging (kg) :	0,652
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	60
Humidity min	(N/A)	Humidity max	(N/A)



30/11/2022 21:11 - This article is available on our website : www.gsh-med.fr with the following reference HDACCE128N