

CASE DORO FOR 2414/2424, BLACK



Réf. HDACCE130



CASE DORO FOR 2414/2424, BLACK



Réf. HDACCE130

FEATURES			
Medical device	No		
Marquage		Barcode	7322460075255
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	8,0
Length (cm) :		Length packaging (cm) :	15,5
Height (cm) :		Height packaging (cm) :	3,0
Weight (kg) :		Weight packaging (kg) :	0,040
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	60
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:07 - This article is available on our website : www.gsh-med.fr with the following reference HDACCE130