

FLIP CASE FOR DORO 8080



Réf. HDACCE135



FLIP CASE FOR DORO 8080



Réf. HDACCE135

FEATURES			
Medical device	No		
Marquage		Barcode	7322460076528
Guarantee (month):	0	Customs Code	4202390090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	12,5
Length (cm) :		Length packaging (cm) :	17,5
Height (cm) :		Height packaging (cm) :	2,0
Weight (kg) :		Weight packaging (kg) :	0,174
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	60
Humidity min	(N/A)	Humidity max	(N/A)



19/10/2022 12:38 - This article is available on our website : www.gsh-med.fr with the following reference HDACCE135