

OXYGEN EMPTY BACKPACK



Réf. 580825



OXYGEN EMPTY BACKPACK



Réf. 580825

Bags and cases for transporting emergency response equipment. Can be used to safely carry medical oxygen cylinders and accessories. Designed for all cylinder sizes offered by Spengler.

FEATURES			
Medical device	No		
Marquage		Barcode	3700446034641
Guarantee (month):	0	Customs Code	90200000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	3,500
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



15/06/2022 21:05 - This article is available on our website : www.gsh-med.fr with the following reference 580825

All characteristics are given for information only and may be changed