

PRESSURE REDUCER WITH BRACKET 0-15L/MIN

Réf. 840123



PRESSURE REDUCER WITH BRACKET

0-15L/MIN



Réf. 840123

Gas release devices. Medical gas regulators reduce the pressure of medical gases and are used, when equipped with a flowmeter, to regulate the flow of medical gas administered to patients. They are connected to cylinders containing high-pressure medical gases, equipped with a tap. Product technological lifespan : 10 years. Full compulsory product inspection by a certified technician within 5 years, 3 years if under intensive use. Regular product inspection recommended yearly.

FEATURES			
Medical device	Yes	Classe IIb	
Marquage	CE2460	Barcode	8592346325241
Guarantee (month):	12	Customs Code	90192000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,680
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



15/06/2022 21:06 - This article is available on our website : www.gsh-med.fr with the following reference 840123

All characteristics are given for information only and may be changed