

## VINYL RESUSCITATORS ADULT (MASK 5)



Réf. 810102



CE  
0120



## VINYL RESUSCITATORS ADULT (MASK 5)



Réf. 810102

Emergencies - resuscitation - first aid. Flexible, lightweight and ergonomic, these resuscitators maintain safe and effective ventilation, as much in cases requiring resuscitation as in cases where breathing is low. These systems can be connected to an oxygen supply.

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0120	Barcode	5030267073276
Guarantee (month):	0	Customs Code	90192000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,590
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



15/06/2022 21:06 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 810102

All characteristics are given for information only and may be changed