

## VINYL MASK SINGLE USE, ADULT, YELLOW, T3

Réf. 860132



CE  
0120



# VINYL MASK SINGLE USE, ADULT, YELLOW, T3



**Réf. 860132**

Emergencies - resuscitation - first aid. Can be fitted to the manual insufflators. These facial resuscitation masks ensure administration of oxygen and anesthetic vapours from anesthesia administration devices. They can be used in spontaneously breathing patients or to ensure controlled ventilation. Their specially designed shape means they fit all morphologies.

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0120	Barcode	5030267079698
Guarantee (month):	0	Customs Code	90192000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,030
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



15/06/2022 21:07 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 860132

All characteristics are given for information only and may be changed