


PRINTER FOR TOURNIQUET 800 SERIE



Réf. GHP552885



PRINTER FOR TOURNIQUET 800 SERIE



Réf. **GHP552885**

Accessory for electronic tourniquets

FEATURES			
Medical device	No		
Marquage	Barcode	8715946465784	
Guarantee (month):	0	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	2,530
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:05 - This article is available on our website : www.gsh-med.fr with the following reference GHP552885

All characteristics are given for information only and may be changed