


LABEL ROLL FOR PRINTER SERIE 800

Hpm™

Réf. GHP552886



LABEL ROLL FOR PRINTER SERIE 800



Réf. **GHP552886**

Accessory for electronic tourniquets

FEATURES			
Medical device	No		
Marquage		Barcode	3700446034276
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,344
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:06 - This article is available on our website : www.gsh-med.fr with the following reference GHP552886