

DISPOSABLE OXYGEN MASK WITH TUBING HIGH CONCENTRATION PEDIATRIC



Réf. 580873

CE
0120



DISPOSABLE OXYGEN MASK WITH TUBING HIGH CONCENTRATION PEDIATRIC



Réf. **580873**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0120	Barcode	5030267051779
Guarantee (month):	0	Customs Code	90192000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	15,0
Length (cm) :		Length packaging (cm) :	15,0
Height (cm) :		Height packaging (cm) :	6,0
Weight (kg) :		Weight packaging (kg) :	0,068
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



15/06/2022 21:06 - This article is available on our website : www.gsh-med.fr with the following reference 580873

All characteristics are given for information only and may be changed