

DISPOSABLE SINGLE CUFF 61X10CM YELLOW (CARTON OF 10 PIECES)

Réf. GHP552894

HpmTM



CE
0123



DISPOSABLE SINGLE CUFF 61X10CM YELLOW (CARTON OF 10 PIECES)



Réf. **GHP552894**

ETO Sterilized. 2 m long hose made of TPU. Straight and Conical. No Latex. Wide range of sizes. Color coded. Customized protection sleeve included. CPC male as standard. Soft and easy to apply

FEATURES			
Medical device	Yes	Classe Is (fonction stérile)	Accuracy measurement
Marquage	CE0123		Classe I
Guarantee (month):	12	Barcode	7350054600231
Sterile	No	Customs Code	
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	22,0
Length (cm) :		Length packaging (cm) :	55,0
Height (cm) :		Height packaging (cm) :	21,0
Weight (kg) :		Weight packaging (kg) :	2,390
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:06 - This article is available on our website : www.gsh-med.fr with the following reference GHP552894

All characteristics are given for information only and may be changed