

DISPOSABLE SINGLE CUFF 109X13CM GREEN (CARTON OF 10 PIECES)

Réf. GHP552898

Hpm™



CE
0123



Disposable Tourniquet Cuff
PN 15-180 Size 110x10cm 43"
DiattoSolut2 S.p.a.
Scandinavian Medical Systems AB
Västberga Årå 368
SE-126 50 Hågersten
Sweden

STERILE EO 2013-04-10 LOT 1104001

Welling Ltd 016 Welling Tech Consulting
Suite B, 29 Harley Street London W1G 9QR
England, United Kingdom
Tel: +44(20) 88184372, Fax: +44(20)76811847

CE 0123 Reg. 002010-12-20

DISPOSABLE SINGLE CUFF 109X13CM GREEN (CARTON OF 10 PIECES)



Réf. **GHP552898**

ETO Sterilized. 2 m long hose made of TPU. Straight and Conical. No Latex. Wide range of sizes. Color coded. Customized protection sleeve included. CPC male as standard. Soft and easy to apply

FEATURES			
Medical device	Yes	Classe Is (fonction stérile)	Accuracy measurement Classe I
Marquage	CE0123		Barcode 7350054600279
Guarantee (month):	12		Customs Code
Sterile	No		
INFORMATIONS		Main unit	Packaged unit
Width (cm) :		Width packaging (cm):	26,0
Length (cm) :		Length packaging (cm) :	64,0
Height (cm) :		Height packaging (cm) :	20,5
Weight (kg) :		Weight packaging (kg) :	3,858
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:06 - This article is available on our website : www.gsh-med.fr with the following reference GHP552898