

DISPOSABLE SINGLE CONICAL CUFF 61X10CM RED (CARTON OF 10 PIECES)

Réf. GHP552900

Hpm™



CE
0123



Disposable Tourniquet Cuff
PIN 75-140 Size 63x9cm, 25"
Distributed by:
Scandinavian Medical Systems AB
Västberga Ålls 36B
SE-120 50 Hågersten
Sweden
STERILE/10
2013-04-10
LOT
1104001
Welling Ltd (a Welling Tech Consulting
Suite B, 29 Manly Street London W10 3QR
England, United Kingdom
Tel: +44(0)20 88168312, Fax: +44(0)20 758 11843
CE

DISPOSABLE SINGLE CONICAL CUFF 61X10CM RED (CARTON OF 10 PIECES)



Réf. **GHP552900**

ETO Sterilized. 2 m long hose made of TPU. Straight and Conical. No Latex. Wide range of sizes. Color coded. Customized protection sleeve included. CPC male as standard. Soft and easy to apply

FEATURES			
Medical device	Yes	Classe Is (fonction stérile)	Accuracy measurement Classe I
Marquage	CE0123		Barcode 7350054600350
Guarantee (month):	12		Customs Code
Sterile	No		
INFORMATIONS		Main unit	Packaged unit
Width (cm) :		Width packaging (cm) :	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) : 2,650	
SPECIFIC STORAGE CONDITION			
Temperature min (N/A)		Temperature max (N/A)	
Humidity min (N/A)		Humidity max (N/A)	



31/03/2022 02:06 - This article is available on our website : www.gsh-med.fr with the following reference GHP552900

All characteristics are given for information only and may be changed