

DISPOSABLE SINGLE CONICAL CUFF 76X10CM YELLOW (CARTON OF 10 PIECES)

Réf. GHP552901

Hpm™



CE
0123



Disposable Tourniquet Cuff
PIN 75-150 Size 75x10cm 30"
Distributed by:
Scandinavian Medical Systems AB
Västbergsgatan 36B
SE-120 50 Hägersten
Sweden
STERILE EO LOT 2013-06-10 1104021
CE 0123
Welling Ltd Via Welling Tech Consulting
Suite B, 25 Harley Street London W1G 9DP
England, United Kingdom
Tel: +44(20) 88168372, Fax: +44(20) 76811843

DISPOSABLE SINGLE CONICAL CUFF 76X10CM YELLOW (CARTON OF 10 PIECES)



Réf. **GHP552901**

ETO Sterilized. 2 m long hose made of TPU. Straight and Conical. No Latex. Wide range of sizes. Color coded. Customized protection sleeve included. CPC male as standard. Soft and easy to apply

FEATURES			
Medical device	Yes	Classe Is (fonction stérile)	Accuracy measurement
Marquage	CE0123		Classe I
Guarantee (month):	12	Barcode	7350054600347
Sterile	No	Customs Code	
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	27,0
Length (cm) :		Length packaging (cm) :	64,0
Height (cm) :		Height packaging (cm) :	20,0
Weight (kg) :		Weight packaging (kg) :	3,174
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:06 - This article is available on our website : www.gsh-med.fr with the following reference GHP552901

All characteristics are given for information only and may be changed