

AUTOCLAVABLE SILICONE CUFF, 52 X 7CM (20.8 INCH), SINGLE

Réf. 554075



AUTOCLAVABLE SILICONE CUFF,52 X 7CM (20.8 INCH), SINGLE



Réf. **554075**

| FEATURES | | | |
|----------------------------|-----------|-------------------------|----------------|
| Medical device | No | | |
| Marquage | CE | Barcode | 28809315676601 |
| Guarantee (month): | 12 | Customs Code | 9018 90 2090 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 7,0 | Width packaging (cm): | |
| Length (cm) : | 52,0 | Length packaging (cm) : | |
| Height (cm) : | | Height packaging (cm) : | |
| Weight (kg) : | | Weight packaging (kg) : | 0,200 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



31/03/2022 02:50 - This article is available on our website : www.gsh-med.fr with the following reference 554075

All characteristics are given for information only and may be changed