

DTS-3000 STAND, 102CM X 54CM

Réf. 554115



DTS-3000 STAND, 102CM X 54CM



Réf. 554115

FEATURES			
Medical device	No		
Marquage	CE	Barcode	18809315676182
Guarantee (month):	12	Customs Code	9018 90 2090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	56,0
Length (cm) :		Length packaging (cm) :	56,0
Height (cm) :		Height packaging (cm) :	23,0
Weight (kg) :		Weight packaging (kg) :	11,600
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:51 - This article is available on our website : www.gsh-med.fr with the following reference 554115

All characteristics are given for information only and may be changed