



ISOTHERM POUCH WITH REUSABLE COLD GEL (WITH VELCRO) DE BOISSY

De Boissy
— DEPUIS 1982

Réf. MDPISOT01



ISOTHERM POUCH WITH REUSABLE COLD GEL (WITH VELCRO) DE BOISSY



Réf. **MDPISOT01**

FEATURES			
Medical device	No		
Marquage	AUCUN	Barcode	3661040219783
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	14,0	Width packaging (cm):	
Length (cm) :	19,0	Length packaging (cm) :	
Height (cm) :	4,0	Height packaging (cm) :	
Weight (kg) :	0,010	Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



23/06/2022 12:31 - This article is available on our website : www.gsh-med.fr with the following reference MDPISOT01

All characteristics are given for information only and may be changed