


## PACK DORO SMARTPHONE 8100 PLUS AND PINK WATCH



Réf. HDSMWAT01RC



## PACK DORO SMARTPHONE 8100 PLUS AND PINK WATCH



Réf. HDSMWAT01RC

FEATURES			
Medical device	No		
Marquage	Barcode		
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



11/11/2022 12:38 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference HDSMWAT01RC

## PACK DORO SMARTPHONE 8100 PLUS AND PINK WATCH



Réf. HDSMWAT01RC

### COMPOSITION

Code	Descriptive	Quantity
HDSMAR01	Smartphone DORO 8100 PLUS, graphite color	1
HDWATC01RC	DORO Watch Pink	1