

CUSHING RETRACTOR 25CM 30MM

Réf. IE08830



CUSHING RETRACTOR 25CM 30MM



Réf. IE08830

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040016443
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	2,7	Width packaging (cm):	
Length (cm) :	24,6	Length packaging (cm) :	
Height (cm) :	2,0	Height packaging (cm) :	
Weight (kg) :	0,059	Weight packaging (kg) :	0,059
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:31 - This article is available on our website : www.gsh-med.fr with the following reference IE08830

All characteristics are given for information only and may be changed