

BANKHART SHOULDER RETRACTOR 24 CM - 50X10 MM

Réf. IE08853



BANKHART SHOULDER RETRACTOR 24 CM - 50X10 MM



Réf. IE08853

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | Yes | Classe IIa | |
| Marquage | CE | Barcode | 3661040084671 |
| Guarantee (month): | 24 | Customs Code | 9018.9090 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 2,5 | Width packaging (cm): | |
| Length (cm) : | 23,8 | Length packaging (cm) : | |
| Height (cm) : | 2,0 | Height packaging (cm) : | |
| Weight (kg) : | 0,108 | Weight packaging (kg) : | 0,108 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



30/03/2022 16:31 - This article is available on our website : www.gsh-med.fr with the following reference IE08853

All characteristics are given for information only and may be changed