

MOHR TUBING CLAMP WITH SCREW 18X18

Réf. IS13618



MOHR TUBING CLAMP WITH SCREW 18X18



Réf. IS13618

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040030777
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,8	Width packaging (cm):	
Length (cm) :	1,8	Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,015	Weight packaging (kg) :	0,015
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:43 - This article is available on our website : www.gsh-med.fr with the following reference IS13618

All characteristics are given for information only and may be changed