

GIGLI SAW 50 CM (SET OF 10)

Réf. IS15100



GIGLI SAW 50 CM (SET OF 10)



Réf. **IS15100**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040030784
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	12,1	Width packaging (cm):	
Length (cm) :	15,3	Length packaging (cm) :	
Height (cm) :	0,6	Height packaging (cm) :	
Weight (kg) :	0,027	Weight packaging (kg) :	0,027
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:43 - This article is available on our website : www.gsh-med.fr with the following reference IS15100

All characteristics are given for information only and may be changed