

CARRYING BAG FOR ASPIRET, ASKIR 20 & ASKIR 30

Réf. BCASPI129

CA-MI



CARRYING BAG FOR ASPIRET, ASKIR 20 & ASKIR 30



Réf. BCASPI129

FEATURES			
Medical device	No		
Marquage		Barcode	3661040165769
Guarantee (month):	12	Customs Code	4202121900
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	21,0	Width packaging (cm):	
Length (cm) :	45,0	Length packaging (cm) :	
Height (cm) :	37,0	Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,650
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	1	Humidity max	85



30/03/2022 23:10 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI129