

BAG ISOTHERMAL SHOULDER STRAP FOR DIABETIQUE ELITE BAGS, BLACK

Réf. MEFIT02_N

 ELITE BAGS

GARANTIE
2
ANS
YEARS
GUARANTEE



BAG ISOTHERMAL SHOULDER STRAP FOR DIABETIQUE ELITE BAGS, BLACK



Réf. MEFIT02_N

FEATURES			
Medical device	No		
Marquage		Barcode	8436564540720
Guarantee (month):	24	Customs Code	4202129990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	14,0	Width packaging (cm):	18,0
Length (cm) :	20,0	Length packaging (cm) :	21,2
Height (cm) :	8,0	Height packaging (cm) :	8,0
Weight (kg) :	0,460	Weight packaging (kg) :	0,432
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:19 - This article is available on our website : www.gsh-med.fr with the following reference MEFIT02_N

All characteristics are given for information only and may be changed