

ENGEL PLASTER SAW 15 CM

Réf. IS15115



ENGEL PLASTER SAW 15 CM



Réf. **IS15115**

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040030814
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	7,2	Width packaging (cm):	7,2
Length (cm) :	15,2	Length packaging (cm) :	15,2
Height (cm) :	0,9	Height packaging (cm) :	0,9
Weight (kg) :	0,105	Weight packaging (kg) :	0,105
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/05/2022 12:35 - This article is available on our website : www.gsh-med.fr with the following reference IS15115

All characteristics are given for information only and may be changed