

HANDLE FOR GIGLI SAW BY PAIR STAINLESS STEEL

Réf. IS15150



HANDLE FOR GIGLI SAW BY PAIR STAINLESS STEEL



Réf. **IS15150**

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040030838
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	7,2	Width packaging (cm):	7,2
Length (cm) :	11,5	Length packaging (cm) :	11,5
Height (cm) :	1,2	Height packaging (cm) :	1,2
Weight (kg) :	0,056	Weight packaging (kg) :	0,056
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:43 - This article is available on our website : www.gsh-med.fr with the following reference IS15150

All characteristics are given for information only and may be changed