

ADSON RETRACTOR 4X5 PRONGS BLUNT 330MM

Réf. IE11235



ADSON RETRACTOR 4X5 PRONGS BLUNT 330MM



Réf. **IE11235**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040044170
Guarantee (month):	24	Customs Code	9018908500
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	10,8	Width packaging (cm):	
Length (cm) :	32,2	Length packaging (cm) :	
Height (cm) :	3,6	Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,184
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:32 - This article is available on our website : www.gsh-med.fr with the following reference IE11235

All characteristics are given for information only and may be changed