

EAR SYRINGE 110 ML

Réf. IS17110



EAR SYRINGE 110 ML



Réf. IS17110

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040030951
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	6,0	Width packaging (cm):	
Length (cm) :	23,0	Length packaging (cm) :	
Height (cm) :	5,5	Height packaging (cm) :	
Weight (kg) :	0,279	Weight packaging (kg) :	0,280
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/05/2022 12:35 - This article is available on our website : www.gsh-med.fr with the following reference IS17110