

## HEADSET FOR MOBILE PHONE 334/338/342/345/410/610/612/680/682

Réf. HDACCE15



# HEADSET FOR MOBILE PHONE

## 334/338/342/345/410/610/612/680/682



Réf. HDACCE15

FEATURES			
Medical device	No		
Marquage		Barcode	7322460051754
Guarantee (month):	24	Customs Code	8518309590
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	7,9	Width packaging (cm):	7,9
Length (cm) :	16,5	Length packaging (cm) :	16,5
Height (cm) :	3,7	Height packaging (cm) :	3,7
Weight (kg) :	0,031	Weight packaging (kg) :	0,031
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	60
Humidity min	(N/A)	Humidity max	(N/A)



30/11/2022 21:10 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference HDACCE15

2/2