

VACHER SNARE

Réf. IS18431



VACHER SNARE



Réf. IS18431

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040087764
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	18,0	Width packaging (cm):	18,0
Length (cm) :	35,0	Length packaging (cm) :	35,0
Height (cm) :	1,0	Height packaging (cm) :	1,0
Weight (kg) :	0,200	Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:44 - This article is available on our website : www.gsh-med.fr with the following reference IS18431

All characteristics are given for information only and may be changed