

HARTMANN EUSTACHIAN CATHETER 3 MM

Réf. IS19103



HARTMANN EUSTACHIAN CATHETER 3 MM



Réf. **IS19103**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040031002
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	7,0	Width packaging (cm):	
Length (cm) :	17,7	Length packaging (cm) :	
Height (cm) :	0,5	Height packaging (cm) :	
Weight (kg) :	0,014	Weight packaging (kg) :	0,014
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:44 - This article is available on our website : www.gsh-med.fr with the following reference IS19103

All characteristics are given for information only and may be changed