

PALMER SPECULUM 5 MM

Réf. IS20105



PALMER SPECULUM 5 MM



Réf. IS20105

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040054643
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,010	Weight packaging (kg) :	0,010
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:45 - This article is available on our website : www.gsh-med.fr with the following reference IS20105

All characteristics are given for information only and may be changed