

SIMS VAGINAL SPECULUM SET OF 3

Réf. IS20300



SIMS VAGINAL SPECULUM SET OF 3

Réf. **IS20300**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040031057
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	9,3	Width packaging (cm):	
Length (cm) :	16,1	Length packaging (cm) :	
Height (cm) :	3,2	Height packaging (cm) :	
Weight (kg) :	0,340	Weight packaging (kg) :	0,340
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:45 - This article is available on our website : www.gsh-med.fr with the following reference IS20300

All characteristics are given for information only and may be changed