


NOGENT TYPE RETRACTOR (WITH 6 PAIRS VALVES) (LASS. P-74)

Réf. IE14250



NOGENT TYPE RETRACTOR (WITH 6 PAIRS VALVES) (LASS. P-74)



Réf. IE14250

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040016740
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	21,5	Width packaging (cm):	
Length (cm) :	36,0	Length packaging (cm) :	
Height (cm) :	10,0	Height packaging (cm) :	
Weight (kg) :	0,837	Weight packaging (kg) :	0,837
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:34 - This article is available on our website : www.gsh-med.fr with the following reference IE14250

All characteristics are given for information only and may be changed