


RICHARDSON RETRACTOR 25 CM 44 X 38 MM

Réf. IE14438



RICHARDSON RETRACTOR 25 CM 44 X 38 MM



Réf. IE14438

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040016764
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	5,0	Width packaging (cm):	
Length (cm) :	25,0	Length packaging (cm) :	
Height (cm) :	3,6	Height packaging (cm) :	
Weight (kg) :	0,099	Weight packaging (kg) :	0,104
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:34 - This article is available on our website : www.gsh-med.fr with the following reference IE14438

All characteristics are given for information only and may be changed