

DEAVER RETRACTOR 50 MM

Réf. IE15450



DEAVER RETRACTOR 50 MM



Réf. IE15450

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040016818
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	10,5	Width packaging (cm):	
Length (cm) :	30,5	Length packaging (cm) :	
Height (cm) :	5,0	Height packaging (cm) :	
Weight (kg) :	0,327	Weight packaging (kg) :	0,327
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:35 - This article is available on our website : www.gsh-med.fr with the following reference IE15450

All characteristics are given for information only and may be changed