

KILLIAN NASAL SPECULUM 75 MM WITH LATERAL SCREW

Réf. IS23075



KILLIAN NASAL SPECULUM 75 MM WITH LATERAL SCREW



Réf. IS23075

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040031163
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	6,5	Width packaging (cm):	
Length (cm) :	16,2	Length packaging (cm) :	
Height (cm) :	8,0	Height packaging (cm) :	
Weight (kg) :	0,100	Weight packaging (kg) :	0,100
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:46 - This article is available on our website : www.gsh-med.fr with the following reference IS23075

All characteristics are given for information only and may be changed