


DEAVER RETRACTOR 75 MM

Réf. IE15475



DEAVER RETRACTOR 75 MM



Réf. IE15475

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040016825
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	12,0	Width packaging (cm):	
Length (cm) :	31,1	Length packaging (cm) :	
Height (cm) :	7,5	Height packaging (cm) :	
Weight (kg) :	0,452	Weight packaging (kg) :	0,452
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:35 - This article is available on our website : www.gsh-med.fr with the following reference IE15475

All characteristics are given for information only and may be changed