

KILLIAN NASAL SPECULUM 30 MM WITHOUT SCREW

Réf. IS23130



KILLIAN NASAL SPECULUM 30 MM WITHOUT SCREW



Réf. **IS23130**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040031187
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	8,5	Width packaging (cm):	
Length (cm) :	23,5	Length packaging (cm) :	
Height (cm) :	1,7	Height packaging (cm) :	
Weight (kg) :	0,102	Weight packaging (kg) :	0,102
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:46 - This article is available on our website : www.gsh-med.fr with the following reference IS23130

All characteristics are given for information only and may be changed