

KILLIAN NASAL SPECULUM 50 MM WITHOUT LATERAL SCREW

Réf. IS23250



KILLIAN NASAL SPECULUM 50 MM WITHOUT LATERAL SCREW



Réf. **IS23250**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040031200
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	6,0	Width packaging (cm):	
Length (cm) :	16,0	Length packaging (cm) :	
Height (cm) :	1,0	Height packaging (cm) :	
Weight (kg) :	0,100	Weight packaging (kg) :	0,108
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:46 - This article is available on our website : www.gsh-med.fr with the following reference IS23250

All characteristics are given for information only and may be changed