

COVER FOR JAR 1000 ML AND 2000 ML

CA-MI

Réf. BCASPI84



COVER FOR JAR 1000 ML AND 2000 ML

CA-MI

Réf. BCASPI84

FEATURES			
Medical device	No		
Marquage		Barcode	3661040164786
Guarantee (month):	12	Customs Code	90330090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,200	Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:33 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI84

2/2